

# FluoriMax<sup>®</sup>

2.5% Sodium Fluoride Varnish

With Hydroxyapatite

**Rx Only**

Shake Well Prior to Application

Contents 8 mL

Storage: Store at room temperature. (23° ± 2° C).

Do not freeze or expose to extreme heat



elevate  
ORAL CARE<sup>®</sup>

**FluoriMax<sup>®</sup> Varnish is a thin liquid that is applied differently than other varnishes.**

**Please read directions thoroughly before using.**

### Directions for Use:

1. A prophylaxis is not required, however, removal of excessive plaque or calculus may be beneficial.
2. For the first use from each bottle, replace the storage cap with the flip-top dropper cap.
3. Immediately prior to application, shake the bottle vigorously for 15 seconds.
4. While holding the bottle upright, flip open the dropper cap.
5. Dispense the desired amount into a dappen dish or mixing well and close the cap. As a guide, use these amounts based on the dentition being treated. Adjust as necessary:

Site specific use	1 drop
Primary dentition only	1 drops
Mixed dentition	2 drops
Adult dentition only	3 drops
6. Leave the teeth moist, but immediately before application have the patient swallow or use suction to remove any pooled saliva. If the dentition is dry, have the patient swish with water, then suction.
7. Retract the cheek on one side of the mouth and start with the bottom arch. Dip the brush into the dappen dish and make *one swipe only*, from the posterior to the midline on the buccal/labial surface of one quadrant. Repeat on the top arch. Repeat the procedure on the other side of the mouth.

For primary dentition, have the patient bite together. Allowing the brush bristles to touch both upper and lower arches simultaneously, make *one swipe only* from posterior to the midline on the buccal/labial surface. Repeat the procedure on the other side of the mouth.

For enhanced drying, have the patient breathe in while keeping the cheek retracted for 3-5 seconds.

**Patient Instructions:**

After treatment, instruct the patient to allow FluoriMax to remain on the teeth for at least 4 hours.

1. Avoid alcohol-containing beverages and mouth rinses for at least 4 hours. A normal diet can resume immediately, since unlike sticky and tacky varnishes the FluoriMax varnish layer forms a tenacious tooth bond.
2. Wait 4 - 6 hours before resuming normal oral hygiene habits including brushing or flossing treated teeth. If possible, wait until the next morning.

**Caution:** U.S. Federal Law restricts this device to sale by or on the order of a dental or medical professional.

**Storage:** Store in a cool, dry place at room temperature. (23° ± 2° C) Do not freeze or expose to extreme heat.

**Questions? Comments?**

**Call 877-866-9113 • [www.elevateoralcare.com](http://www.elevateoralcare.com)**

**Flammable:** The product contains ethyl alcohol. Keep away from open flame. Product safe once applied and varnish has adhered to tooth surface.



**Composition:** 1 mL of FluoriMax contains 25 mg sodium fluoride, equivalent to 11.3 mg fluoride ion in an alcohol solution of natural resins. The FluoriMax formula does not include any gluten, pine, peanut, tree-nut, egg, milk, fish, shellfish, soy or wheat products.

**Indications:** A topical fluoride for the treatment of dentinal or post-operative sensitivity.

**Contraindications:** Ulcerative gingivitis and stomatitis.

**Interactions:** The use of other fluoride preparations (e.g. foams, gels or rinses) should not be administered on the same day of application of FluoriMax. The use of fluoride supplements should be discontinued for several days after application.

**Adverse Reactions:** The following has been observed with some fluoride varnishes: Edematous swelling has been reported only in rare instances, especially following application to extensive surfaces. Dyspnea, although extremely rare, has occurred in asthmatic children. Nausea has been reported following extensive applications.

**Dosage:** The fluoride content in FluoriMax is dosed in a way that neither acute nor chronic side effects are to be expected if applied according to instructions.

**Made in U.S.A. by:**

Elevate Oral Care  
West Palm Beach, FL 33411  
877-866-9113

FluoriMax is a registered trademark of Elevate Oral Care, LLC. ELE639-07212020